

Avian Influenza Virus H7N9 Real Time RT-PCR Kit User Manual

LT029030RR

For use with ABI Prism™ 7000/3000/7500/7900/Step One Plus, iCycler iQ™4/iQ™5; Smart Cycler II, Bio-Rad CFX 96; Rotor Gene™ 6000; Mx3000P/3005P; ML-Opticon2/Chromo4; LightCycler™ 480 Instrument

Life Technologies (India) Pvt. Ltd.
306, Aggarwal City Mall, Opposite MZK Pitampura, Delhi - 110094 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechnologies.com Website: www.lifetechnologies.com



1. Intended Use

Avian influenza virus H7N9 real time RT-PCR kit is used for the detection of gene H7 and gene N9 of avian influenza A subtype H7N9 in human nasal and pharyngeal secretions and bird feces by using real time PCR systems.

2. Principle of Real-Time PCR

The principle of the real-time detection is based on the fluorogenic 5' nuclease assay. During the PCR reaction, the DNA polymerase cleaves the probe at the 5' end and separates the reporter dye from the quencher dye only when the probe hybridizes to the target DNA. This cleavage results in the fluorescent signal generated by the cleaved reporter dye, which is monitored real-time by the PCR detection system. The PCR cycle at which an increase in the fluorescence signal is detected initially (Ct) is proportional to the amount of the specific PCR product. Monitoring the fluorescence intensities during Real Time allows the detection of the accumulating product without having to re-open the reaction tube after the amplification.

3. Product Description

Highly pathogenic avian influenza (HPAI) caused by certain subtypes of influenza A virus in animal populations, particularly chickens, poses a continuing global human public health risk. Direct human infection by an avian influenza A (H5N1) virus was first recognized during the 1997 outbreak in Hong Kong. The avian influenza virus H7N9 is one subgroup among the larger group of H7 viruses. Some cases of human infection with H7N9 virus in China are confirmed till early April of 2013.

Avian influenza virus H7N9 real time RT-PCR kit contains a specific ready-to-use system for the detection of avian influenza virus H7N9 by Reverse Transcription Polymerase Chain Reaction (RT-PCR) in the real-time PCR system. The master contains Super Mix for the specific amplification of the avian influenza virus H7N9 RNA. The reaction is done in one step real time RT-PCR. The first step is a reverse transcription (RT), during which the avian influenza virus H9 RNA is transcribed into cDNA. Afterwards, a thermostable DNA polymerase is used to amplify the specific gene fragments by polymerase chain reaction. Fluorescence is emitted and measured by the real time systems' optical unit during the PCR. The detection of amplified avian influenza virus H7N9 DNA fragment is performed in fluorimeter channel FAM and HEX/VIC/JOE with the fluorescent quencher BHQ1. In addition, the kit contains a system to identify possible PCR inhibition by measuring the Cal Red 610/ROX/TEXAS RED fluorescence of the internal control (IC).

4. Kit Contents

Ref.	Type of reagent	Presentation	25rxns
1	H7N9 Super Mix	1 vial, 480µl	
2	RT-PCR Enzyme Mix	1 vial, 28µl	
3	Molecular Grade Water	1 vial, 400µl	
4	H7N9 Internal Control	1 vial, 30µl	
5	H7N9 Positive Control	1 vial, 30µl	

Analysis sensitivity: 1 X 10³ copies/ml;

Note: Analysis sensitivity depends on the sample volume, elution volume, nucleic acid extraction methods and other factors. If you use the RNA extraction kits recommended, the analysis sensitivity is the same as it declares. However, when the sample volume is dozens or even hundreds of times greater than elution volume by some concentrating method, it can be much.

5. Storage

- All reagents should be stored at -20°C. Storage at +4°C is not recommended.
- All reagents can be used until the expiration date indicated on the kit label.
- Repeated thawing and freezing (> 3x) should be avoided, as this may reduce the sensitivity of the assay.
- Cool all reagents during the working steps.
- Super Mix should be stored in the dark.

6. Additionally Required Materials and Devices

- Biological cabinet
- Vortex mixer
- Cryo-container
- Sterile filter tips for micro pipets
- Disposable gloves, powderless
- Refrigerator and Freezer
- Desktop microcentrifuge for "ependorf" type tubes (RCF max. 16,000 x g)
- Real time PCR system
- Real time PCR reaction tubes/plates
- Pipets (0.5µl – 1000µl)
- Sterile microtubes
- Biohazard waste container
- Tube racks

7. Warnings and Precaution

- Carefully read this instruction before starting the procedure.
- For in vitro diagnostic use only.
- This assay needs to be carried out by skilled personnel.
- Clinical samples should be regarded as potentially infectious materials and should be prepared in a laminar flow hood.
- This assay needs to be run according to Good Laboratory Practice.
- Do not use the kit after its expiration date.
- Avoid repeated thawing and freezing of the reagents, this may reduce the sensitivity of the test.
- Once the reagents have been thawed, vortex and centrifuge briefly the tubes before use.
- Prepare quickly the Reaction mix on ice or in the cooling block.
- Set up two separate working areas: 1) Isolation of the RNA/ DNA and 2) Amplification/ detection of amplification products.
- Pipets, vials and other working materials should not circulate among working units.
- Use always sterile pipette tips with filters.
- Wear separate coats and gloves in each area.
- Do not pipette by mouth. Do not eat, drink, smoke in laboratory.
- Avoid aerosols

8. Sample Collection, Storage and transport

- Collected samples in sterile tubes;
- Specimens can be extracted immediately or frozen at -20°C to -80°C.
- Transportation of clinical specimens must comply with local regulations for the transport of etiologic agents

9. Procedure

9.1 RNA-Extraction

Different brand RNA extraction kits are available. You may use your own extraction systems or the commercial kit based on the yield. For the RNA extraction, please comply with the manufacturer's instructions. The recommended Extraction kit is as follows:


RNA Extraction Kit GEN 52-904 LT

9.2 Internal Control

It is necessary to add internal control (IC) in the reaction mix. Internal Control (IC) allows the user to determine and control the possibility of PCR inhibition. Add the internal control (IC) 1µl/rxn and the result will be shown in the Cal Red 610/ROX/TEXAS RED.

9.3 RT-PCR Protocol

The Master Mix volume for each reaction should be pipetted as follows:



※PCR system without Cal Red 610/ROX/TEXAS RED channel may be treated with 1µl Molecular Grade Water instead of 1µl IC.

- The volumes of Super Mix and Enzyme Mix per reaction multiply with the number of samples, which includes the number of controls, standards, and sample prepared. Molecular Grade Water is used as the negative control. For reasons of unprecise pipetting, always add an extra virtual sample. Mix completely then spin down briefly in a centrifuge.
- Pipet 20µl Master Mix with micropipets of sterile filter tips to each of the Real time PCR reaction plate/tubes. Separately add 5µl RNA sample, positive and negative controls to different reaction plate/tubes. Immediately close the plate/tubes to avoid contamination.
- Spin down briefly in order to collect the Master Mix in the bottom of the reaction tubes.
- Perform the following protocol in the instrument:

Temperature/Time	Cycle	Selection of fluorescence channels
45°C for 10min	1cycle	FAM
95°C for 15min	1cycle	HEX/VIC/JOE
95°C for 15sec, 60°C for 1min	45cycles	Cal Red 610/ROX/TEXAS RED

- If you use ABI Prism® system, please choose "none" as passive reference and quencher.

10. Threshold setting: just above the maximum level of molecular grade water.

11. Quality control: Negative control, internal control and positive control must be performed correctly, otherwise the sample results are invalid.

Control	Channel			Ct value	
	FAM	HEX/VIC/JOE	Cal Red 610/ROX/TEXAS RED		
Molecular Grade Water	UNDET	UNDET	UNDET	25-35	
Positive Control	≤35	≤35	—		

12. Data Analysis and Interpretation : The following results are possible:

	Ct value			Result Analysis
	FAM	HEX	Cal Red 610	
1#	UNDET	UNDET	25-35	Below the detection limit or negative
2#	≤43	UNDET	—	Gene H7 positive;
3#	UNDET	≤43	—	Gene N9 positive;
4#	≤43	≤43	—	H7N9 Positive;
5#	43~45	—	25-35	Re-test; If it is still 43~45, report as 1#
6#	UNDET	UNDET	UNDET	PCR Inhibition; No diagnosis can be concluded.