

## *n*Taq-Tenuto (*UDG plus*)

Cat.#	Size	Size
P225AU	250 units	5 units/ $\mu$ l
P250AU	500 units	5 units/ $\mu$ l

Store at -20°C

Supplied with: 10X *n*Taq-Tenuto Buffer (Mg<sup>2+</sup> plus)  
dNTP Mixture with dUTP (2 mM each)  
Sterile water

### India Contact:

**Life Technologies (India) Pvt. Ltd.**  
306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).  
Ph: +91-11-42208000, 42208111, 42208222  
Mobile: +91-9810521400  
Fax: +91-11-42208444  
Email: [customerservice@lifetechnindia.com](mailto:customerservice@lifetechnindia.com)  
Web: [www.lifetechnindia.com](http://www.lifetechnindia.com)

### Product description

*n*Taq-Tenuto (*UDG plus*) is *n*Taq (Cat.# P025, P050) supplemented with 3'→5' proofreading activity and a PCR enhancing factor for improved efficiency and fidelity. *n*Taq-Tenuto (*UDG plus*) can be used to amplify DNA longer than 10 kb, which is difficult with common Taq polymerases alone. Thus, this product is improved in both fidelity (> 2 fold) of PCR products and amplification efficiency of longer PCR products. Also, *UDG* and dUTP are included in the mixture to prevent the reamplification of carryover PCR products between reactions.

### Characteristics

- Carry-over contamination control: contains *UDG*
- Molecular weight: 94 kDa
- Error rate: 3.0 X 10<sup>-6</sup>
- Thermal stability: Half life of 40 min at 95°C
- A-tail formation at 3' ends of amplified DNA products.

### Applications

- Amplification of long DNA fragments (>5~15 kb)
  - Amplification of high-complexity template DNA such as cDNA and genomic DNA
  - Primer extension
  - Colony PCR
  - Multiplex PCR
  - Labeling of DNA fragments with radioactive-isotopes
  - Nucleotide sequencing
- For Research Use Only. Not for use in diagnostic procedures.*

ISO9001 ISO14001 ISO13485

## *n*Taq-Tenuto (*UDG plus*)

Cat.#	Size	Size
P225AU	250 units	5 units/ $\mu$ l
P250AU	500 units	5 units/ $\mu$ l

Store at -20°C

Supplied with: 10X *n*Taq-Tenuto Buffer (Mg<sup>2+</sup> plus)  
dNTP Mixture with dUTP (2 mM each)  
Sterile water

### India Contact:

**Life Technologies (India) Pvt. Ltd.**  
306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).  
Ph: +91-11-42208000, 42208111, 42208222  
Mobile: +91-9810521400  
Fax: +91-11-42208444  
Email: [customerservice@lifetechnindia.com](mailto:customerservice@lifetechnindia.com)  
Web: [www.lifetechnindia.com](http://www.lifetechnindia.com)

*For Research Use Only. Not for use in diagnostic procedures.*

ISO9001 ISO14001 ISO13485

### Product description

*n*Taq-Tenuto (*UDG plus*) is *n*Taq (Cat.# P025, P050) supplemented with 3'→5' proofreading activity and a PCR enhancing factor for improved efficiency and fidelity. *n*Taq-Tenuto (*UDG plus*) can be used to amplify DNA longer than 10 kb, which is difficult with common Taq polymerases alone. Thus, this product is improved in both fidelity (> 2 fold) of PCR products and amplification efficiency of longer PCR products. Also, *UDG* and dUTP are included in the mixture to prevent the reamplification of carryover PCR products between reactions.

### Characteristics

- Carry-over contamination control: contains *UDG*
- Molecular weight: 94 kDa
- Error rate: 3.0 X 10<sup>-6</sup>
- Thermal stability: Half life of 40 min at 95°C
- A-tail formation at 3' ends of amplified DNA products.

### Applications

- Amplification of long DNA fragments (>5~15 kb)
- Amplification of high-complexity template DNA such as cDNA and genomic DNA
- Primer extension
- Colony PCR
- Multiplex PCR
- Labeling of DNA fragments with radioactive-isotopes
- Nucleotide sequencing

### Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase free
- Inhibitor-free

### Unit definition

One unit is defined as the amount of enzyme required to incorporate 10 nmol of dNTP into acid insoluble materials in 30 min at 74°C in a 50- $\mu$ l reaction mixture (20 mM Tris-HCl/pH 8.8, 50 mM KCl, 2.5 mM MgCl<sub>2</sub>, 10 mM  $\beta$ -mercaptoethanol, 12.5  $\mu$ g of calf thymus DNA).

### Storage buffer

20 mM Tris-HCl (pH 7.9), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% NP-40, 0.5% Tween-20, 50% glycerol.

### Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase free
- Inhibitor-free

### Unit definition

One unit is defined as the amount of enzyme required to incorporate 10 nmol of dNTP into acid insoluble materials in 30 min at 74°C in a 50- $\mu$ l reaction mixture (20 mM Tris-HCl/pH 8.8, 50 mM KCl, 2.5 mM MgCl<sub>2</sub>, 10 mM  $\beta$ -mercaptoethanol, 12.5  $\mu$ g of calf thymus DNA).

### Storage buffer

20 mM Tris-HCl (pH 7.9), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% NP-40, 0.5% Tween-20, 50% glycerol.

Your Molecular & Cell Technology Partner

### Standard PCR conditions

- PCR mixture <sup>a</sup>	
10X <i>n</i> Taq-Tenuto Buffer (Mg <sup>2+</sup> plus)	2 $\mu$ l
<i>n</i> Taq-Tenuto ( <i>UDG plus</i> ) <sup>b</sup> (5 units/ $\mu$ l)	0.2 $\mu$ l
dNTP mixture with dUTP(2 mM each, final conc., 200 $\mu$ M each)	2 $\mu$ l
Template DNA <sup>c</sup> (0.1~500 ng/ $\mu$ l)	1 $\mu$ l
Primer 1 (5 pmole/ $\mu$ l)	1 $\mu$ l
Primer 2 (5 pmole/ $\mu$ l)	1 $\mu$ l
Sterile water	up to 20 $\mu$ l
<sup>a</sup> Assemble the reaction mixture on ice	
<sup>b</sup> Add the PCR polymerase at the final step	
<sup>c</sup> Plasmid DNA, 0.1 ng~30 ng; genomic DNA, 50 ng~500 ng	

### - PCR cycle

Pre-incubation (for <i>UDG</i> )	25°C	10 min
Initial denaturation	95°C	2 min
Denaturation	95°C	30 sec
Annealing <sup>a</sup>	55°C~65°C	30~60 sec
Elongation	72°C	1 min/kb
Number of cycles	25~35 times	
Final elongation	72°C	5 min

When cycles are over, keep the reaction mixture at 4°C; may add 10 mM EDTA until use to prevent DNA degradation.  
<sup>a</sup>Recommended annealing temperatures is 5 to 10°C below the lower T<sub>m</sub> of the two primers used.

Your Molecular & Cell Technology Partner

### Standard PCR conditions

- PCR mixture <sup>a</sup>	
10X <i>n</i> Taq-Tenuto Buffer (Mg <sup>2+</sup> plus)	2 $\mu$ l
<i>n</i> Taq-Tenuto ( <i>UDG plus</i> ) <sup>b</sup> (5 units/ $\mu$ l)	0.2 $\mu$ l
dNTP mixture with dUTP(2 mM each, final conc., 200 $\mu$ M each)	2 $\mu$ l
Template DNA <sup>c</sup> (0.1~500 ng/ $\mu$ l)	1 $\mu$ l
Primer 1 (5 pmole/ $\mu$ l)	1 $\mu$ l
Primer 2 (5 pmole/ $\mu$ l)	1 $\mu$ l
Sterile water	up to 20 $\mu$ l
<sup>a</sup> Assemble the reaction mixture on ice	
<sup>b</sup> Add the PCR polymerase at the final step	
<sup>c</sup> Plasmid DNA, 0.1 ng~30 ng; genomic DNA, 50 ng~500 ng	

### - PCR cycle

Pre-incubation (for <i>UDG</i> )	25°C	10 min
Initial denaturation	95°C	2 min
Denaturation	95°C	30 sec
Annealing <sup>a</sup>	55°C~65°C	30~60 sec
Elongation	72°C	1 min/kb
Number of cycles	25~35 times	
Final elongation	72°C	5 min

When cycles are over, keep the reaction mixture at 4°C; may add 10 mM EDTA until use to prevent DNA degradation.  
<sup>a</sup>Recommended annealing temperatures is 5 to 10°C below the lower T<sub>m</sub> of the two primers used.