

BrightBox Quantitation Assay

The BrightBox Assay is a novel method to quantitate your NGS libraries

Library quantitation can be a lengthy, costly, and inaccurate process where variability can lead to unequal library loading and waste time and money.

The BrightBox Assay is a quick and simple method to accurately quantitate any NGS library.

Get results quickly by shortening your current library preparation and instrument loading workflows.



Why Use BrightBox

✓ FAST

Save time with reactions done in just 5 minutes.

✓ EASY

Reduce manual errors and save time with premade Assay Mix and prediluted standards kept in refrigerated storage.

✓ ACCURATE

Improve data with consistent pooling based on actual molarity.

✓ FLEXIBLE

Quantitate any library with attached P5 and P7 adapters.

Faster than current workflows.
Preparation and incubation times for the BrightBox assay take a fraction of the time required for other methods.

- Aliquot Assay Mix
- Add Libraries
- Incubate and read - 5 minutes!

BrightBox

15 minutes

26 libraries per 96 well-plate

- Mix Library with Dye
- Load one sample at a time onto Qubit

Mass Analysis

25 minutes

16 libraries

- Estimate Mass - 1 hour
- Dilute Library
 - Prime
- Run on instrument - 1 hour

Fragment Analysis

2 hours

11 - 96 libraries

- Dilute Libraries - 1 hour
- Thaw Master Mix
- Assemble Master Mix
- Aliquot Master Mix
 - Add Libraries
- Cycle and read - 1 hour



qPCR

2.5 hours

12 libraries per 96 well-plate


Quantitation is correlated to read numbers.

Molarity was calculated for 13 libraries with the BrightBox assay and a commonly used qPCR library quantitation kit. Equal volumes of each library were loaded onto a MiSeq and read numbers were plotted against calculated molarity.




Reproducible results.

16 different libraries were repeatedly run using the BrightBox assay on a qPCR instrument.




Quantitation is linear and reproducible.

Standards from the BrightBox kit were run 12 different times on a qPCR instrument and show high linearity and low variability.



Correlation between methods.

16 different libraries were quantitated with the BrightBox assay and a qPCR library quantitation kit. Outlying points are due to inaccurate estimation of "average" fragment size during qPCR quantitation.



For research use only. Not for use in diagnostic procedures.