

India Contact:

Life Technologies (India) Pvt. Ltd.  
Mobile: +91-9810521400, Ph: +91-11-42208000  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)  
Web: [www.lifetechindia.com](http://www.lifetechindia.com)

## Product Information

### Human Peripheral Blood CD4+ Helper T Cells – Untouched by negative selection

Catalog Number	10HU-023N	Cell Number	1.0 x 10 <sup>7</sup> cells/vial
Species	<i>Homo sapiens</i>	Storage Temperature	Liquid Nitrogen

## Description

The CD4+ T cells, also known as T helper cells (T<sub>h</sub> cells), are a type of T cells that play an important role in the immune system, particularly in the adaptive immune system [1]. CD4+ T cells can regulate immune response through different cytokines secreted [2]. They are essential in B cell antibody class switching, in the activation and growth of cytotoxic T cells, and in maximizing bactericidal activity of phagocytes such as macrophages.

Considering the diverse and important role CD4+ T cells play in the immune system, it is not surprising that these cells often influence the immune response against disease. CD4+ T cells have been involved in hypersensitivity and as the main target cells of HIV pathogenesis.

**iXCells Biotechnologies** offers CD4+ T Cells isolated from normal human peripheral blood mononuclear cells (PBMCs) using negative immunomagnetic selection. > 90% of the cells are CD4+ as showed by flow cytometric analysis.

## Product Details

<b>Tissue</b>	Normal human peripheral blood
<b>Package Size</b>	1.0x10 <sup>7</sup> cells/vial
<b>Passage Number</b>	P0
<b>Shipped</b>	Cryopreserved
<b>Storage</b>	Liquid nitrogen
<b>Growth Properties</b>	Suspension
<b>Media</b>	Blood Cell Culture Medium (Cat# MD-0007)

## Protocols

### Thawing of Frozen Cells

1. Upon receipt of the frozen cells, it is recommended to thaw the cells and initiate the culture immediately in order to retain the highest cell viability.
2. To thaw the cells, put the vial in 37°C water bath with gentle agitation for 1-2 minutes. Keep the cap out of water to minimize the risk of contamination.
3. Pipette the cells into a 15 mL conical tube with 5 mL fresh Blood Cell Culture Medium (Cat# MD-0007).
4. Centrifuge at 400-450 g for 5 minutes under room temperature.
5. Remove the supernatant and cell is ready for downstream applications.

**Safety Precaution:** *it is highly recommended that protective gloves and clothing should be used when handling frozen vials.*

## References

[1] Hu, J.; Paul, W. E. (2008). "CD4 T cells: Fates, functions, and faults". *Blood* **112** (5): 1557.

[2] Toscano MA, Bianco GA, Ilarregui JM, Croci DO, Correale J, Hernandez JD, Zwirner NW, Poirier F, Riley EM, Baum LG, Rabinovich GA. (2007). "Differential glycosylation of TH1, TH2 and TH-17 effector cells selectively regulates susceptibility to cell death". *Nat Immunol* (8): 825–34.

## Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans. While iXCells Biotechnologies uses reasonable efforts to include accurate and up-to-date information on this product sheet, we make no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. iXCells Biotechnologies does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. iXCells Biotechnologies is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, iXCells Biotechnologies is not liable for damages arising from the misidentification or misrepresentation of cultures.

© iXCells Biotechnologies 2015. All rights reserved.