



India Contact:

Life Technologies (India) Pvt. Ltd. Mobile: +91-9810521400, Ph: +91-11-42208000 Email: <u>customerservice@lifetechindia.com</u> Web: <u>www.lifetechindia.com</u>

Product Information

Mouse Hepatic Stellate Cells (MHSC)

Catalog Number	10MU-021	Cell Number	0.5 million cells/vial
Species	Mus Musculus	Storage Temperature	Liquid nitrogen

Description

Hepatic stellate cells (HSCs) are liver-specific mesenchymal cells, and account for 5~8% of the cells in the liver. HSCs play vital roles in the homeostasis of liver extracellular matrix, repair, regeneration and fibrosis, and control retinol metabolism, storage and release. The stellate cell is the major cell type involved in liver fibrosis in response to liver injury ^[1]. In healthy liver, HSCs are in a quiescent state, and contains numerous vitamin A lipid droplets, constituting the largest reservoir of vitamin A in the body. When the liver is damaged, HSCs can change into an activated state, which is characterized by proliferation, contractility and chemotaxis. The amount of vitamin A decreases progressively in injured liver. The activated HSCs also secrete collagen scar tissue, which can lead to cirrhosis. In chronic liver disease, prolonged and repeated activation of stellate cells causes liver fibrosis ^[2]. Primary culture of HSC is a valuable tool to study liver fibrosis.

iXCells Biotechnologies provides high quality frozen or freshly plated mouse HSCs (MHSC) isolated from mouse liver. The cells are frozen immediately after isolation, with >0.5 million cells per vial. There will be 30-50% dead cells after thawing, which is normal for MHSC. We don't recommend passaging the MHSCs.



Figure legend: (A) Phase contrast image of day 3 culture of Mouse Hepatic Stellate Cells (MHSC). (B) Vitamin A droplets in MHSC were visualized by blue autofluorescence (lower panel) when excited with UV light.

1

All Rights Reserved

Product Details

Tissue	Mouse liver
Package Size	0.5 million cells/vial (enough to seed on 3 wells of a 12-well plate).
Passage Number	P0
Shipped	Fresh or frozen
Storage	Liquid nitrogen
Growth Properties	Adherent
Media	Stellate Cell Growth Medium (Cat# MD-0014)

Protocols

- Add warm Stellate Cell Growth Medium (Cat# MD-0014) to culture dish, then thaw frozen cells quickly in water bath, add 1 mL warm culture medium, mix cells gently by pipetting up and down 5 times (avoid vigorous pipetting), then add cells directly to culture dish containing the medium.
- Recommended seeding density: 2~3 million is enough to plate for one 96 well plate, or one 6 well plate, or one 12 well plate.
- **3.** To enhance cell attachment, change medium 24~30 hours after the initial plating. Pipet off the medium gently using pipette to avoid cell detachment during medium change.
- 4. Prior to treatment, starve the cells overnight without growth supplement. Cells can be used for treatment 30 hours after the initial plating.

Note: We don't recommend passaging the MHSCs.

References

[1] Chunyue Yin, et al, (2013) "Hepatic stellate cells in liver development, regeneration, and cancer". J Clin Invest. 123(5):1902–1910.

[2] Rockey D. C. (2001) "Hepatic blood flow regulation by stellate cells in normal and injured liver". Semin Liver Dis 21(3):337-49.

Disclaimers

2

This product is intended for laboratory research purposes only. It is not intended for use in humans. While iXCells Biotechnologies uses reasonable efforts to include accurate and up-to-date information on this product sheet, we make no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. iXCells Biotechnologies does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. iXCells Biotechnologies is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, iXCells Biotechnologies is not liable for damages arising from the misidentification or misrepresentation of cultures. © iXCells Biotechnologies 2015. All rights reserved.

For Research Only

All Rights Reserved

iXCells Biotechnologies USA, Inc.