

## Coating Antibody for Tomato ringspot virus (ToRSV)



Test Format: Double Antibody Sandwich (DAS) ELISA

## Intended Use:

This capture antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

## Includes:

- Capture Antibody
- 96-well microtiter plates, solid

## Option

Coating antibody for Tomato ringspot virus (ToRSV), 500 testwells CAB 22001/0500

Coating antibody for Tomato ringspot virus (ToRSV), 1000 testwells CAB 22001/1000

Coating antibody for Tomato ringspot virus (ToRSV), 5000 testwells CAB 22001/5000