



Coating Antibody for Broad bean wilt virus 1,2 (BBWV-1,2)



Click to view a larger image

Test Format: Double Antibody Sandwich (DAS) ELISA

Intended Use:

This capture antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

Includes:

- · Capture Antibody
- · 96-well microtiter plates, solid

Option

Coating antibody for Broad bean wilt virus 1,2 (BBWV-1,2), 500 testwells CAB 46202/0500

Coating antibody for Broad bean wilt virus 1,2 (BBWV-1,2), 1000 testwells CAB 46202/1000

Coating antibody for Broad bean wilt virus 1,2 (BBWV-1,2), 5000 testwells CAB 46202/5000

India Contact: