



Coating Antibody for Erwinia amylovora (Ea)

Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This capture antibody is intended for use with its corresponding detection antibodies (ECA).

Test performance using other detection antibodies has not been evaluated and is not guaranteed.

Includes:

- Capture Antibody
- 96-well microtiter plates, solid

Option

Coating Antibody for Erwinia Amylovora (Ea), 500 testwells CAB 52101/0500

Coating Antibody for Erwinia Amylovora (Ea), 1000 testwells CAB 52101/1000

Coating Antibody for Erwinia Amylovora (Ea), 5000 testwells CAB 52101/5000