

## Coating Antibody for Potato virus A (PVA)



Click to view a larger image

Test Format: Triple Antibody Sandwich (TAS) ELISA

## Intended Use:

This capture antibody is intended for use with its corresponding detection antibodies (ECA).

Test performance using other detection antibodies has not been evaluated and is not guaranteed.

## Includes:

- Capture Antibody
- 96-well microtiter plates, solid

## Option

Coating antibody for Potato virus A (PVA), 500 testwells CAB 60000/0500

Coating antibody for Potato virus A (PVA), 1000 testwells CAB 60000/1000

Coating antibody for Potato virus A (PVA), 5000 testwells CAB 60000/5000