





Coating Antibody for Pectobacterium atrosepticum (Patro)

Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This capture antibody is intended for use with its corresponding detection antibodies (ECA).

Test performance using other detection antibodies has not been evaluated and is not guaranteed.

Includes:

- Capture Antibody
- 96-well microtiter plates, solid

Option

Coating antibody for Pectobacterium atrosepticum (Patro), 500 testwells CAB 80003/0500

Coating antibody for Pectobacterium atrosepticum (Patro), 1000 testwells CAB 80003/1000

Coating antibody for Pectobacterium atrosepticum (Patro), 5000 testwells CAB 80003/5000

India Contact: