

Coating Antibody for Blueberry leaf mottle virus (BLMoV)



Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This capture antibody is intended for use with its corresponding detection antibodies (ECA).

Test performance using other detection antibodies has not been evaluated and is not guaranteed.

Includes:

- Capture Antibody
- 96-well microtiter plates, solid

Option

Coating antibody for Blueberry leaf mottle virus (BLMoV), 500 testwells CAB 88001/0500

Coating antibody for Blueberry leaf mottle virus (BLMoV), 1000 testwells CAB 88001/1000

Coating antibody for Blueberry leaf mottle virus (BLMoV), 5000 testwells CAB 88001/5000