

Coating Antibody for Pepper mild mottle virus (PMMoV)



Test Format: Double Antibody Sandwich (DAS) ELISA

Intended Use:

This capture antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

Includes:

- Capture Antibody
- 96-well microtiter plates, solid

	Option
	Agdia Coating Antibody PMMoV, 500 testwells CAB 37502/0500
<input type="checkbox"/>	Agdia Coating Antibody PMMoV, 1000 testwells CAB 37502/1000
	Agdia Coating Antibody PMMoV, 5000 testwells CAB 37502/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com