

Enzyme Conjugate for Tobacco streak virus (TSV)



Test Label: Alkaline Phosphatase
Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This detection antibody (bottle A) and enzyme conjugate (bottle B) are intended for use with their corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- **User Guide**

<input type="checkbox"/>	Option
<input type="checkbox"/>	Enzyme conjugate for Tobacco streak virus (TSV), alkphos, 500 testwells ECA 25500/0500
<input type="checkbox"/>	Enzyme conjugate for Tobacco streak virus (TSV), alkphos, 1000 testwells ECA 25500/1000
<input type="checkbox"/>	Enzyme conjugate for Tobacco streak virus (TSV), alkphos, 5000 testwells ECA 25500/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com