

Enzyme Conjugate for Blueberry shoestring virus (BSSV)

Test Label: Horseradish Peroxidase

Test Format: Double Antibody Sandwich (DAS) ELISA

Intended Use:

This enzyme conjugate is intended for use with its corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Peroxidase enzyme conjugate
- MRS component, 5X
- TMB peroxidase substrate solution, 1X
- [User Guide](#)

Option

Enzyme conjugate for Blueberry shoestring virus (BSSV), perox, 500 testwells
ECP 19000/0500

Enzyme conjugate for Blueberry shoestring virus (BSSV), perox, 1000 testwells
ECP 19000/1000

Enzyme conjugate for Blueberry shoestring virus (BSSV), perox, 5000 testwells
ECP 19000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com