

ELISA Reagent Set for Pepper mottle virus (PepMoV)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

The PepMoV ELISA is an effective tool for the detection of *Pepper mottle virus* in plant foliage.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Pepper mottle virus (PepMoV), 96 testwells, SRA 38000/0096

Reagent Set for Pepper mottle virus (PepMoV), 500 testwells, SRA 38000/0500

Reagent Set for Pepper mottle virus (PepMoV), 1000 testwells, SRA 38000/1000

Reagent Set for Pepper mottle virus (PepMoV), 5000 testwells, SRA 38000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com