



## ELISA Reagent Set for CMV II

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal Detection Reagent: Monoclonal

## Intended Use:

This test is a qualitative serological assay for detection of *Cucumber mosaic virus*-subgroup II in plant foliage.

This test system is specific for CMV subgroup II (CMV-ToRS). It does not detect CMV subgroup I (CMV-DTL) or other cucumoviruses including *Peanut stunt virus* and *Tomato aspermy virus*.

Some well-known, characterized isolates that are detected by this test are: Subgroup II (CMV-ToRS): LS, R, R76B, Rm, S, U.

## **Includes:**

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- User Guide
- Tech Note





Reagent Set for Cucumber mosaic virus subgroup II (CMV II), 96 testwells,SRA 44800/0096

Reagent Set for Cucumber mosaic virus subgroup II (CMV II), 500 testwells, SRA 44800/0500

Reagent Set for Cucumber mosaic virus subgroup II (CMV II), 1000 testwells, SRA 44800/1000

Reagent Set for Cucumber mosaic virus subgroup II (CMV II), 5000 testwells, SRA 44800/5000