


Certificate of Analysis

For Laboratory Use Only

PRODUCT NAME	Ginsenoside Rb2
PART NUMBER	00007196
STANDARD TYPE	Analytical Standard (AS)
LOT NUMBER	00007196-301
REPORT NUMBER	CDXA-RSS-4934-02
SAMPLE NUMBER	CDXA-12-2055
DATE OF SAMPLE	04/04/2012
DATE OF RE-EVALUATION	07/01/2016 (1 st); 04/13/2020 (2 nd)
DATE OF REPORT	06/01/2020



CHEMICAL NAMES	(3 β ,12 β)-20-[[6-O- α -L-Arabinopyranosyl- β -D-glucopyranosyl]oxy]-12-hydroxydammar-24-en-3-yl 2-O- β -D-glucopyranosyl- β -D-glucopyranoside; Ginsenoside C
CHEMICAL FORMULA	C ₅₃ H ₉₀ O ₂₂
MOLECULAR WEIGHT (MW)	1079.27
CHEMICAL FAMILY	Saponins
CAS NUMBER	[11021-13-9]
EC#(EINECS)	234-251-4
RTECS	LZ5779240

ANALYTICAL RESULTS

TEST	METHOD	SPECIFICATION	RESULT
HPLC	0.700.10.2.METH102	NA	96.2%
Mass Spectrum	Direct Infusion, ESI (+)	Conforms	Conforms
Appearance	NA	NA	White powder

STORAGE CONDITIONS

STORAGE	Room Temperature in a dry place.
EXPIRATION DATE	04/2025 under the above conditions.

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.


Certificate of Analysis

For Laboratory Use Only

ANALYTICAL CONDITIONS

INSTRUMENT	AGILENT 1260 HPLC UV-VIS (DAD) CONNECTED TO DIONEX CORONA VEO CHARGED AEROSOL DETECTOR (CAD); AGILENT 6510 QTOF
COLUMN	Phenomenex Kinetex 2.6 µm particle size; 150 x 4.6 mm
MOBILE PHASE	A - Ultrapure water, B - Acetonitrile; Start at 15% B then increasing to 30% B over 15 minutes, then increasing to 32% B over one minute, then increasing to 38% B over 3 minutes, then increasing to 43% B over 5 minutes, then increasing to 55% B over 3 minutes, then increasing to 65% B over 4 minutes, then increasing to 70% B over 4 minutes and then increasing to 90% B over 3 minutes. Decrease to 15% B over one minute and then hold for 4 minutes.
COLUMN TEMP.	40 °C
FLOW RATE	1.0 mL/minute
INJECTION VOL.	5.0 µL
INJECTION CONC.	1.2 mg/mL in Methanol
DETECTION	CAD


HPLC CHROMATOGRAM OF GINSENOSE RB2 (CDXA-12-2055)




This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

Certificate of Analysis For Laboratory Use Only

HPLC CHROMATOGRAM OF DILUENT BLANK



MASS SPECTRUM OF GINSENOSE RB2 (CDXA-12-2055)



This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

Certificate of Analysis

For Laboratory Use Only

REVISION HISTORY

<u>Revision History</u>	<u>Date of Revision</u>	<u>Document/Changes</u>
00	04/23/2012	New report
01	07/01/2016	Updated expiration date, updated standard type
02	06/01/2020	Passed second re-evaluation by HPLC. Updated the structure, chemical names, expiration date, analytical results, analytical conditions, and chromatograms. Removed the melting point.

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.