





Certificate of Analysis

For Laboratory Use Only

PRODUCT NAME Ursolic acid

PART NUMBER 00021600

STANDARD TYPE Primary (P)

LOT NUMBER 00021600-022

REPORT NUMBER CDXA-RSS-3985-03


SAMPLE NUMBER CDXA-10-2823

DATE OF SAMPLE 05/05/2010

DATE OF RE-EVALUATION 06/09/2014 (1st); 08/15/2018 (2nd);

07/07/2022 (3rd)

DATE OF REPORT 08/22/2022



CHEMICAL NAMES (3β)-3-Hydroxyurs-12-en-28-oic acid; Bungeolic acid; Formosolic acid; Malol;

Prunol; Urson

CHEMICAL FORMULA C₃₀H₄₈O₃

MOLECULAR WEIGHT (MW) 456.70

CHEMICAL FAMILY Triterpenes

CAS NUMBER [77-52-1]

EC#(EINECS) 201-034-0

RTECS YU9520000

ANALYTICAL RESULTS

TEST	METHOD	SPECIFICATION	RESULT
Adjusted Purity	NA	≥ 85%	95.1%
HPLC	0.700.10.2.METH22	≥ 90%	96.7%
NMR	0.700.12.4	Conforms to structure	Conforms to structure
Mass Spectrum	Direct Infusion, ESI (–)	Conforms	Conforms
Residual Solvent*	99.1-CD-3.0-000115	NA	EtOH -1.2%
Water	0.700.12.37	NA	0.5%
Appearance	NA NA White po		White powder

*Testing performed at an ISO 17025 accredited subcontracted laboratory.

ADJUSTED PURITY: 95.1% IS BASED ON (100 - 1.2 SOLVENTS - 0.5 WATER) X 96.7% HPLC**STORAGE CONDITIONS**

STORAGE Room Temperature in a dry place.

EXPIRATION DATE 07/2027 under the above conditions.

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com







Certificate of Analysis

For Laboratory Use Only

ANALYTICAL CONDITIONS

INSTRUMENT AGILENT 1260 HPLC UV-VIS (DAD) CONNECTED TO DIONEX CORONA VEO

CHARGED AEROSOL DETECTOR (CAD); AGILENT 1100 ION TRAP

COLUMN Phenomenex Luna C18(2) 150 x 2.00 mm, 5 μm particle size

MOBILE PHASE A - 0.1% Trifluoroacetic acid in Ultrapure water, B - Acetonitrile;

Isocratic 5% B for 5 minutes, then increasing to 95% B over 20 minutes. Hold at

95% B for 15 minutes, then decreasing to 5% B over 1 minute.

COLUMN TEMP. 60 °C


FLOW RATE 0.4 mL/minute

INJECTION VOL. 2.0 μL

INJECTION CONC. 1.0 mg/mL in Methanol

DETECTION CAD

HPLC CHROMATOGRAM OF URSOLIC ACID (CDXA-10-2823)











Certificate of Analysis

For Laboratory Use Only

HPLC CHROMATOGRAM OF DILUENT BLANK



MASS SPECTRUM OF URSOLIC ACID (CDXA-10-2823)











Certificate of Analysis

For Laboratory Use Only

¹H-NMR SPECTRUM OF URSOLIC ACID (CDXA-10-2823) IN PYRIDINE-d₅



¹³C-NMR SPECTRUM OF URSOLIC ACID (CDXA-10-2823) IN PYRIDINE-d₅










Certificate of Analysis

For Laboratory Use Only

COSY-NMR SPECTRUM OF URSOLIC ACID (CDXA-10-2823) IN PYRIDINE-d5



REVISION HISTORY

<u>Revision History</u>	<u>Date of Revision</u>	Document/Changes
00	05/20/2010	New report
01	06/12/2014	Passed reevaluation; updated expiration date
02	09/22/2018	Passed second reevaluation. Updated analytical results, analytical conditions, chemical structure format, expiration date and chromatograms. Removed melting point.
03	08/22/2022	Passed third re-evaluation by HPLC, KF, NMR. Updated expiration date, analytical results, analytical conditions, and chromatograms.