





For Laboratory Use Only

SPECIMEN

Certificate of Analysis

PRODUCT Urtica dioica (Stinging nettle) Leaf

PART NUMBER 00030449

REFERENCE TYPE Biomass Reference Material (BRM) *

LOT NUMBER 00030449-108

COMMON NAME Stinging Nettle

Urtica dioica L. [Urticaceae] **LATIN NAME**

PLANT PART Leaf

CDXP NUMBER CDXA-BRM-108

REPORT NUMBER CDXA-BRMR-029-03

01/27/2004 **DATE OF SAMPLE**

DATE OF RE-EVALUATION 01/30/2012 (1st); 04/25/2016 (2nd); 04/20/2020 (3rd)

DATE OF REPORT 06/10/2020

ANALYTICAL RESULTS

TEST	METHOD	RESULT
Appearance	Macroscopy	(1) Milled: Olive green powder
Anatomy/Morphology	Microscopy	(1) Large cysoliths
HP-TLC	CDXA-TLCM-049-01	Conforms

STORAGE CONDITIONS

STORAGE 20-30 °C; Dry storage area; Insect free; Volatile free

EXPIRATION DATE 04/2025 under the above conditions

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

^{*}Note - Biomass Reference Material (BRM) is not a voucher specimen








SPECIMEN


For Laboratory Use Only

Certificate of Analysis

MACROSCOPY



MICROSCOPY



This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.







SPECIMEN

Certificate of Analysis

For Laboratory Use Only

HP-TLC CONDITIONS

STATIONARY PHASE Silica gel 60, F₂₅₄, 10 x 10 cm HP-TLC plates

SAMPLE PREPARATION ~0.5 g + 5 mL methanol, sonicated in hot water bath for ~10 minutes.

MOBILE PHASE Ethyl acetate/ formic acid/ acetic acid/ Milli-Q water


[10/1.1/1.1/2.6]


CHAMBER TEMPERATURE Ambient


DETECTION (1) UV 366 nm

(2) Natural Products reagent → UV 366 nm
 (3) Natural Products reagent → visible light

HP-TLC PLATES







YELLOW LINE = SAMPLE ORIGIN
RED LINE = SOLVENT FRONT @ 70 mm

HP-TLC LANE APPLICATIONS

Lane	ID	Lane	ID
1	MeOH blank 2.0 µL	5	Urtica dioica 6.0 µL
2	Urtica dioica 1.0 µL	6	Urtica dioica 8.0 µL
3	Urtica dioica 2.0 µL	7	Urtica dioica 10.0 µL
4	Urtica dioica 4.0 µL	8	Urtica dioica 12.0 µL

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.







SPECIMEN

For Laboratory Use Only

Certificate of Analysis

REVISION HISTORY

<u>Revision History</u>	<u>Date of Revision</u>	<u>Document/Changes</u>
00	06/18/2007	New report
01	01/30/2012	Passed re-evaluation; updated expiration date; general template update; updated TLC conditions; updated TLC lane applications; updated TLC method; updated TLC pictures
02	04/25/2016	Passed re-evaluation; updated expiration date; general template update; updated TLC conditions; updated TLC lane applications; updated TLC pictures; updated TLC pictures
03	06/10/2020	Material passed reevaluation by HPTLC. Updated template, expiration date, TLC conditions, plate images, and lane assignments.

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.