

HUMAN UMBILICAL VEIN ENDOTHELIAL CELLS DISSOCIATION SOLUTION

Catalog Number: hUVEC-DM-10

Product description:

Dissociation solution is necessary for gentle detachment of cells from the growth surface. This solution does not contain any protein. Growth media should be added to stop the dissociation activity, further cells can be used for subculturing.

Specifications:

Quantity	10mL
Application	Detaching of Human umbilical vein endothelial cells
Product use	This product is for research use only.
Storage	Store in -20 °C
Shelf life	3 months at -20 °C
Shipping	Dry ice
QC	Tested negative bacteria, fungi, yeast. Tested negative for mycoplasma, endotoxin.

Instructions to use:

Thaw dissociation solution at 37 °C before use.

Note: Dissociation solution once opened, should be used within a month.

Life Technologies (India) Pvt Ltd.

Life Technologies® is a registered trademark of Life Technologies (India) Pvt Ltd.