

Cat No. L-100T1KB

Product Description

Life Technologies India DNA Ladder L-100T1KB is a ready to use mix of 11 proprietary DNA fragments ranging from 100 to 1500 base pairs, suitable for use as molecular standards for agarose gel electrophoresis. The mix contains Orange G and Xylene Cyanol FF for tracking. The mix is recommended to be stored at 4°C for use within 12 months and should be stored at -20°C for use within 24 months.

Concentration: 100 µg per ml

Package: 50 µg per 500ul

Recommended Load: 5µl per well

