

Cat No. L-100T3KB

Product Description

Life Technologies India DNA Ladder L-100T3KB is a ready to use mix of 12 proprietary DNA fragments ranging from 100 to 3000 base pairs, suitable for use as molecular standards for agarose gel electrophoresis. The mix contains Orange G and Xylene Cyanol FF for tracking. The mix is recommended to be stored at 4°C for use within 12 months and should be stored at -20°C for use within 24 months.

Concentration: 108 µg per ml

Package: 54 µg per 500ul

Recommended Load: 5µl per well

