

Product Description

iXCells Biotechnologies provides high quality Human Iliac Artery Endothelial Cells (HIAEC), which are isolated from human iliac artery. HIAEC are shipped in proliferating culture (P3) with >90% confluence in T25 flask, which contains >5 x 10⁵ cells. HIAEC express vWF/Factor VIII and CD31 (PECAM), and DiI-Ac-LDL by uptake. They are negative for HIV-1, HBV, HCV, mycoplasma, bacteria, yeast and fungi and can further expand for 12 population doublings in Endothelial Cell Growth Medium (Cat# MD-0010) under the condition suggested by iXCells Biotechnologies.

Product Details

Tissue	Human iliac artery
Package Size	>90% confluence in T25 flask (5 x 10 ⁵ cells)
Passage Number	P3
Shipped	Proliferating culture
Storage	N/A
Growth Properties	Adherent
Media	Endothelial Cell Growth Medium (Cat# MD-0010)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com