

Product Description

iXCells Biotechnologies provides high quality Human Internal Thoracic Artery Endothelial Cells (HITAEC), which are isolated from normal human internal thoracic (aka mammary) artery and cryopreserved at P2, with >0.5 million cells in each vial. HITAEC express vWF/Factor VIII and CD31 (PECAM), and Dil-Ac-LDL by uptake. They are negative for HIV-1, HBV, HCV, mycoplasma, bacteria, yeast, and fungi and can further expand for 10 population doublings in Endothelial Cell Growth Medium (Cat# MD-0010) under the condition suggested by iXCells Biotechnologies.

Product Details

Tissue	Human internal thoracic (aka mammary) artery
Package Size	5 x 10 ⁵ cells/vial
Passage Number	P2
Shipped	Proliferating culture
Storage	Cryopreserved
Growth Properties	Adherent
Media	Endothelial Cell Growth Medium (Cat# MD-0010)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com