

Product Description

The keratocytes, or corneal fibroblasts, are highly specialized cells that are sandwiched between orthogonally arranged layers of collagen lamellae in the corneal stroma. They play a key role in maintaining the structure and transparency of the cornea, as they are the source of stromal collagen and proglycans. They also play important roles in corneal wound healing and tissue repair, and are known to undergo phenotypic transformations in wounds due to the influence of growth factors and cytokines [1]. Under normal conditions, the keratocytes in the adult cornea are relatively quiescent cells. In the event of corneal injury or trauma, however, the keratocytes differentiate into active, synthesizing cells and rapidly replace damaged stromal matrix. Cultured human keratocytes express functional IL-4Rs [2] and IL-17R [3] on the cell surface, suggesting that these cells may contribute to the role of IL-4 and IL-17 as mediators of allergic reactions in the cornea. Changes in gene expression were observed in keratocytes after interleukin-1 treatment, which provides important insight into gene expression and suggests novel therapeutic targets for the control of corneal inflammation.

iXCells Biotechnologies provides high quality Human Keratocytes (HK), which are isolated from human cornea and cryopreserved at P0, with >0.5 million cells in each vial. HK express fibronectin and are characterized by fibroblast morphology. They are negative for HIV-1, HBV, HCV, mycoplasma, bacteria, yeast, and fungi and can further expand for 12 population doublings in Fibroblast Growth Medium (Cat# MD-0011) under the condition suggested by iXCells Biotechnologies.

Product Details

Tissue	Human cornea
Package Size	0.5 million cells/vial
Passage Number	P0
Shipped	Cryopreserved
Storage	Liquid nitrogen
Growth Properties	Adherent
Media	Fibroblast Growth Medium (Cat# MD-0011)

References

- [1] Fini, M. E. (1999) Keratocyte and fibroblast phenotype in the repairing cornea. *Prog Retin Eye Res* 18:529-551.
- [2] Fukuda, K., Fujitsu, Y., Kumagai, N., Nishida, T. (2002) Characterization of the interleukin-4 receptor complex in human corneal fibroblasts. *Invest Ophthalmol. Vis. Sci.* 43(1):183-8.
- [3] Maertzdorf, J., Osterhaus, A. D., Verjans, G. M. (2002) IL-17 expression in human herpetic stromal keratitis: modulatory effects on chemokine production by corneal fibroblasts. *J. Immunol.* 169(10):5897-903.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
 Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com