


Human TGF-β1 ELISA Kit

PRODUCT SPECIFICATIONS

Assay Range	31.25 pg/ml - 2000 pg/ml
Sensitivity	9 pg/ml
Target species	Human
Incubation	4h30
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 µl
Cross Reaction	Cross reactivity with sheep and pig TGF-β1
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure with Safety Data Sheet
Synonym	TGF-b1, TGF-beta1



TGF-b1 ELISA Kit

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com