


Human VEGF-A ELISA Kit

PRODUCT SPECIFICATIONS

Assay Range	31.25 pg/ml - 1000 pg/ml
Sensitivity	7.9 pg/ml
Target species	Human
Specificity	Specific of human Vascular Endothelial Growth Factor A
Incubation	4h30
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	50 µl
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .
Synonym	VPF



VEGF-A ELISA Kit

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com