


Human p53 ELISA Kit

PRODUCT SPECIFICATIONS

Assay Range	3.12 U/ml - 100 U/ml
Sensitivity	1.5 U/ml
Target species	Human
Specificity	Recognizes both natural and recombinant soluble human p53
Incubation	3h45
Sample Type	Serum, Plasma, Cell lysate
Sample Size	100 µl
Cross Reaction	No cross reactivity with other human soluble molecules
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .



Example of Human p53 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com