


Human CD116 ELISA kit

PRODUCT SPECIFICATIONS

Assay Range	31.25 pg/ml - 1000 pg/ml
Sensitivity	16 pg/ml
Target species	Human
Specificity	Recognizes both natural and recombinant human soluble GM-CSF receptor a chain
Incubation	3h45
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 µl
Cross Reaction	No cross reactivity with other human soluble receptors
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .
Synonym	GM-CSF R



Example of Human CD116 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com