


Human CD178 ELISA kit

PRODUCT SPECIFICATIONS

Assay Range	62.5 pg/ml - 2000 pg/ml
Sensitivity	12 pg/ml
Target species	Human
Specificity	Recognizes both natural and recombinant soluble human Fas ligand (FasL)
Incubation	3h45
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 µl
Cross Reaction	No cross reactivity with other human soluble molecules
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .
Synonym	FAS-Ligand, CD95L



Example of Human CD178 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com