


Human CD14 ELISA kit

PRODUCT SPECIFICATIONS

Assay Range	6.87 ng/ml - 220 ng/ml
Sensitivity	6 ng/ml
Target species	Human
Specificity	Recognizes both natural and recombinant human soluble CD14
Incubation	3h45
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 µl
Cross Reaction	No cross reactivity with other human soluble molecules
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .



Example of Human CD14 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com