


# Human Perforin ELISA Kit

## PRODUCT SPECIFICATIONS

Assay Range	62.5 pg/ml - 2000 pg/ml
Sensitivity	40 pg/ml
Target species	Human
Specificity	Recognizes natural human Perforin
Incubation	2h45
Sample Type	Supernatant, Buffered Solutions, Not suitable for use with serum and plasma.
Sample Size	100 µl
Cross Reaction	No cross reactivity with Trail, TRAIL R3, TRAIL R4, Granzyme B, TNFα, IFNγ, Fas, Fas L, IL-1b, IL-2
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a <a href="#">detailed procedure</a> .



Example of Human Perforin Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)