

# Human IL-10 ELISA Set

## PRODUCT SPECIFICATIONS

Assay Range	12.5 pg/ml - 400 pg/ml
Sensitivity	5 pg/ml
Target species	Human
Incubation	From sample to end 2h45
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100µl
Cross Reaction	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
Kit Content	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, including buffer composition. <a href="#">detailed procedure</a> and <a href="#">Safety Data Sheet</a>

### India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)