

# Human IL-1 $\beta$ ELISA Set

## PRODUCT SPECIFICATIONS

Assay Range	15.6 pg/ml - 500 pg/ml
Sensitivity	5 pg/ml
Target species	Human
Specificity	Recognizes both natural (pro and mature) and recombinant human IL-1b
Incubation	From sample to end 3h45
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 $\mu$ l
Cross Reaction	No cross reactivity with other human cytokines
Kit Content	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and <a href="#">detailed procedure</a> including buffer composition.
Synonym	IL-1b, IL-1beta

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)