

Human IL-17A ELISA Set

PRODUCT SPECIFICATIONS

| | |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Assay Range | 3.12 pg/ml - 100 pg/ml |
| | |
| Sensitivity | 3 pg/ml |
| | |
| Target species | Human |
| | |
| Incubation | From sample to end 3h45 |
| | |
| Sample Type | Serum, Plasma, Cell culture supernatant |
| | |
| Sample Size | 100µl |
| | |
| Cross Reaction | No cross reactivity with other human cytokines |
| | |
| Kit Content | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, including buffer composition. detailed procedure and Safety Data Sheet |
| | |
| | |

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com