


# Human IL-15 ELISA Kit

## PRODUCT SPECIFICATIONS

|                |  |
|----------------|--|
| Assay Range    | 31.25 pg/ml - 1000 pg/ml   |
| Sensitivity    | 18.3 pg/ml   |
| Target species | Human  |
| Specificity    | Recognizes both natural and recombinant human IL-15  |
| Incubation     | 3h45   |
| Sample Type    | Serum, Plasma, Cell culture supernatant  |
| Sample Size    | 50 µl  |
| Cross Reaction | No cross reaction with IL-2, IL-4, IL-7, IL-9, IL-10, IL-13, IL-21 and IFN-g   |
| Kit Content    | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a <a href="#">detailed procedure</a> . |



Example of Human IL-15 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)