


Human IL-33 ELISA Kit

PRODUCT SPECIFICATIONS

Assay Range	31.25 pg/ml - 1000 pg/ml
Sensitivity	12.2 pg/ml
Target species	Human
Specificity	Recognizes both natural and recombinant human IL-33
Incubation	3h45
Sample Type	Serum, Plasma and other body fluids, Cell culture supernatants, Buffered solutions
Sample Size	100 µl (50µl for serum/plasma)
Cross Reaction	No cross reactivity with IL-1b, IL-2, IL-1a, IL-10, IL-12p70, Gp130, IL-17A, IL-23, IFNg, TNFa
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .



Example of Human IL-33 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com