


# Human MCP-1 ELISA kit

## PRODUCT SPECIFICATIONS

Assay Range	31.25 pg/ml - 1000 pg/ml
Sensitivity	5.8pg/ml
Target species	Human
Specificity	Recognizes both natural and recombinant human MCP-1
Incubation	2h15
Sample Type	Serum, plasma, cell culture supernatant, buffered solutions and other body fluids.
Sample Size	100 µl (50µl for serum/plasma)
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a <a href="#">detailed procedure</a> .
Synonym	CCL2, SCYA-2, JE



Example of Human MCP-1 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)