


# Human Galectin-9 ELISA Kit

## PRODUCT SPECIFICATIONS

Assay Range	15.6-500 pg/ml
Sensitivity	8 pg/ml
Target species	Human
Incubation	3h45
Sample Type	Cell culture supernatant, serum, plasma
Sample Size	50µl / 100µl
Cross Reaction	No cross reactivity with GM-CSF, MICB, IL-34, IFNg, IL-29, Galectin-1, Galectin-3, Galectin-8, Galectin-10
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a <a href="#">detailed procedure</a> .
Synonym	HUAT, Ecalectin, LGALS9



Example of Human Galectin-9 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)