


# Human IL-6 ELISA Kit

## PRODUCT SPECIFICATIONS

Assay Range	6.25 pg/ml - 200 pg/ml
Sensitivity	2 pg/ml
Target species	Human
Specificity	Recognizes both natural and recombinant human IL-6
Incubation	1h45
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 µl
Cross Reaction	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a <a href="#">detailed procedure</a> .



Example of Human IL-6 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)