


Human IL-8 ELISA Kit

PRODUCT SPECIFICATIONS

Assay Range	31.25 pg/ml - 1000 pg/ml
Sensitivity	12.3 pg/ml
Target species	Human
Specificity	Recognizes both natural and recombinant human IL-8
Incubation	2h15
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 µl
Cross Reaction	No cross reactivity with other human cytokines
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .
Synonym	CXCL8



Example of Human IL-8 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com