


Human IL-10 ELISA Kit

PRODUCT SPECIFICATIONS

Assay Range	12.5 pg/ml - 400 pg/ml
Sensitivity	4.9 pg/ml
Target species	Human
Specificity	Recognizes both natural and recombinant human IL-10
Incubation	2h45
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 µl
Cross Reaction	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .



Example of Human IL-10 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com