


Human TNF- α ELISA Kit

PRODUCT SPECIFICATIONS

Assay Range	25 pg/ml - 800 pg/ml
Sensitivity	8 pg/ml
Target species	Human
Specificity	Recognizes both natural and recombinant human TNF- α
Incubation	3h45
Sample Type	Serum, Plasma, Cell culture supernatant
Sample Size	100 μ l
Cross Reaction	No cross reactivity with other human cytokines. Cross reactivity with macac TNF- α . No interference with soluble TNFR I and TNFR II.
Kit Content	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure .
Synonym	TNF- α , TNF-alpha



Example of Human TNF- α Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com