


Human CD25 ELISA kit

PRODUCT SPECIFICATIONS

| | |
|----------------|--|
| Assay Range | 68.75 pg/ml - 2200 pg/ml |
| Sensitivity | 32.5 pg/ml |
| Target species | Human |
| Specificity | Recognizes both natural and recombinant human soluble IL-2R |
| Incubation | 3h45 |
| Sample Type | Serum, Plasma, Cell culture supernatant |
| Sample Size | 100 µl |
| Cross Reaction | No cross reactivity with other human soluble receptors |
| Kit Content | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure . |
| Synonym | IL-2R |



Example of Human CD25 Standard Curve

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000 Email: customerservice@lifetechindia.com