

BAPTA, tetrasodium salt


Catalog #: 50004

Product Description

BAPTA and its derivatives are calcium chelators that are commonly used to form calcium buffers with well-defined calcium concentrations. By injecting the chelators into cells or by incubating cells with the AM ester form of the chelators, one can control the cytosolic calcium concentration, an important means to study the roles of calcium. Key advantages of these calcium chelators include relative insensitivity toward intracellular pH change and fast release of calcium. Biotium offers several BAPTA chelators with calcium dissociation constants covering the biologically significant range from 10^{-7} to 10^{-2} M.

- Kd (no Mg^{2+}): 3.6 μ M
- White to off-white solid soluble in water
- Store at 4°C
- $C_{22}H_{18}Br_2K_4N_2O_{10}$
- MW: 786.59

Chemical Structure



Additional Information

Size 100 mg

References

NIL

MSDS (PDF): [MSDS 50004](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, www.lifetechindia.com