

Coating Antibody for Johnsongrass mosaic virus (JGMV)



Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This capture antibody is intended for use with its corresponding detection antibodies (ECA).

Test performance using other detection antibodies has not been evaluated and is not guaranteed.

Includes:

- Capture Antibody
- 96-well microtiter plates, solid

Option

Coating antibody for Johnsongrass mosaic virus (JGMV), 500 testwells
CAB 18300/0500

Coating antibody for Johnsongrass mosaic virus (JGMV), 1000 testwells
CAB 18300/1000

Coating antibody for Johnsongrass mosaic virus (JGMV), 5000 testwells
CAB 18300/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com