

## Detection Antibody for *Xanthomonas axonopodis* pv. *begoniae* (Xab)



[Click to view a larger image](#)

**Test Format:** Antigen coated plate (ACP) ELISA  
**Synonym:** *Xanthomonas campestris* pv. *begoniae*

### Intended Use:

This detection antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

### Includes:

- Detection Antibody
- 96-well microtiter plates, solid

### Option

Detection antibody for *Xanthomonas axonopodis* pv. *begoniae* (Xab), 500 testwells  
CAB 23700/0500

Detection antibody for *Xanthomonas axonopodis* pv. *begoniae* (Xab), 1000 testwells  
CAB 23700/1000

Detection antibody for *Xanthomonas axonopodis* pv. *begoniae* (Xab), 5000 testwells  
CAB 23700/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)