

## Coating Antibody for Potato virus T (PVT)

---



**Test Format:** Triple Antibody Sandwich (TAS) ELISA

### Intended Use:

---

This capture antibody is intended for use with its corresponding detection antibodies (ECA).

Test performance using other detection antibodies has not been evaluated and is not guaranteed.

### Includes:

- Capture Antibody
- 96-well microtiter plates, solid

---

### Option

Coating antibody for Potato virus T (PVT), 500 testwells  
CAB 40100/0500

---

Coating antibody for Potato virus T (PVT), 1000 testwells  
CAB 40100/1000

---

Coating antibody for Potato virus T (PVT), 5000 testwells  
CAB 40100/5000

---

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)