

## Coating Antibody for Groundnut ringspot virus & Tomato chlorotic spot virus (GRSV & TCSV)



**Test Format:** Double Antibody Sandwich (DAS) ELISA

### Intended Use:

This capture antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

### Includes:

- Capture Antibody
- 96-well microtiter plates, solid

### Option

Coating antibody for Groundnut ringspot virus & Tomato chlorotic spot virus (GRSV & TCSV), 500 testwells  
CAB 61000/0500

Coating antibody for Groundnut ringspot virus & Tomato chlorotic spot virus (GRSV & TCSV), 1000 testwells  
CAB 61000/1000

Coating antibody for Groundnut ringspot virus & Tomato chlorotic spot virus (GRSV & TCSV), 5000 testwells  
CAB 61000/5000

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)