

## Coating Antibody for Chrysanthemum virus B (CVB)



**Test Format:** Double Antibody Sandwich (DAS) ELISA

### Intended Use:

This capture antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

### Includes:

- Capture Antibody
- 96-well microtiter plates, solid

### Option

CVB coating antibody, 500 testwells  
CAB 81001/0500

CVB coating antibody, 1000 testwells  
CAB 81001/1000

CVB coating antibody, 5000 testwells  
CAB 81001/5000

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)